

## PPG (Proteose Peptone Glucose Medium)

Freshwater protozoa

<b>Medium</b>	<b>per litre</b>
Proteose peptone (Oxoid L85) *	15.0 g
D-glucose	18.0 g
Page's Amoeba Saline Solution (PAS - see recipe below)	1.0 litre

Mix above constituents and sterilize by pressure cooking at 10 psi for 20 minutes.

### Supply

\* Unipath Ltd, Wade Road, Basingstoke, Hants, RG24 0PW, UK

## PAS (Page's Amoeba Saline Solution)

Constituent of protozoa media

<b>Stocks</b>	<b>per 500 ml</b>
(1) NaCl	12.0 g
MgSO <sub>4</sub> .7H <sub>2</sub> O	0.40 g
CaCl <sub>2</sub> .6H <sub>2</sub> O	0.60 g
(2) Na <sub>2</sub> HPO <sub>4</sub>	14.20 g
KH <sub>2</sub> PO <sub>4</sub>	13.60 g

<b>Medium</b>	<b>per litre</b>
Stock solutions (1) and (2)	5.0 ml each

Add stock solutions to deionized water and make up to 1 litre.