

## Chapman-Andresen's Modified Pringsheim's Solution

Protozoa

<b>Stocks</b>	<b>per 100 ml</b>
(1) Ca(NO <sub>3</sub> ) <sub>2</sub> ·4H <sub>2</sub> O	20.0 g
(2) MgSO <sub>4</sub> ·7H <sub>2</sub> O	2.0 g
(3) Na <sub>2</sub> HPO <sub>4</sub> ·2H <sub>2</sub> O	2.0 g
(4) KCl	2.6 g
(5) FeSO <sub>4</sub> ·7H <sub>2</sub> O	0.2 g
conc H <sub>2</sub> SO <sub>4</sub>	0.1 ml

<b>Medium</b>	<b>per litre</b>
Stock solutions 1 - 5	1.0 ml each

Make up to 1 litre with deionized water. Autoclave at 15 psi for 15 minutes.

Reviewed: 10<sup>th</sup> August 2020