

NN (Non-Nutrient (Amoeba Saline) Agar)

Freshwater protozoa

Stocks (PAS – Page’s Amoeba saline stocks) per 500 ml

(1) NaCl	12.0 g
MgSO ₄ .7H ₂ O	0.40 g
CaCl ₂ .6H ₂ O	0.60 g
(2) Na ₂ HPO ₄	14.20 g
KH ₂ PO ₄	13.60 g

Medium per litre

Stock solutions 1 - 2	5.0 ml each
Bacteriological Agar (Oxoid L11) *	15.0 g

Add stock solutions to deionized water and make up to 1 litre. In suitable vessels thoroughly disperse the agar in cold amoeba saline solution. Autoclave at 15 psi for 15 minutes.

Supply

* Unipath Ltd, Wade Road, Basingstoke, Hants, RG24 0PW, UK

Reviewed: 10th August 2020