

RPL (Rye Grass-Prescott Liquid)

Protozoa

Medium

| | |
|---|------------------|
| Dried milled rye grass * | per litre |
| Prescott's & James's Solution (PJ - see recipe below) | 1.0 g |
| | 1.0 litre |

Bring PJ to the boil then add rye grass. Continue to boil for 5 minutes. Allow to cool then restore volume to 1 litre with deionized water. Filter through GF/C paper and autoclave at 15 psi for 15 minutes.

Supply

* Please enquire

PJ (Prescott's & James's Solution)

Freshwater protozoa

Stocks

| | |
|---|-------------------|
| | per 100 ml |
| (1) $\text{CaCl}_2 \cdot 2\text{H}_2\text{O}$ | 0.43 g |
| KCl | 0.16 g |
| (2) K_2HPO_4 | 0.51 g |
| (3) $\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$ | 0.28 g |

Medium

| | |
|-----------------------|------------------|
| Stock solutions 1 - 3 | per litre |
| | 1.0 ml each |

Make up to 1 litre with deionized water. Autoclave at 15 psi for 15 minutes.