

f/2 + Si (Guillard's medium for diatoms)

Stocks	per 200 ml
(1) NaNO ₃	15 g
(2) NaH ₂ PO ₄ ·2H ₂ O	1.13 g
(3) Trace elements (x10 concentration)	per 200 ml
Na ₂ EDTA	8.32 g
FeCl ₃ ·6H ₂ O	6.30 g
CuSO ₄ ·5H ₂ O	0.02 g
ZnSO ₄ ·7H ₂ O	0.044 g
CoCl ₂ ·6H ₂ O	0.02 g
MnCl ₂ ·4H ₂ O	0.36 g
Na ₂ MoO ₄ ·2H ₂ O	0.012 g
(4) Vitamin mix:	per litre
Cyanocobalamin (Vitamin B ₁₂)	0.0005 g
Biotin	0.0005 g
Thiamine HCl (Vitamin B ₁)	0.1g

For ease of measuring, first prepare primary stock solutions and divide into 1ml aliquots to be frozen for later use.

Component	Primary Stock solution	Quantity	Molar Concentration in final medium
Thiamine HCl (B ₁)	-	100mg	2.96×10^{-7}
Biotin	0.1g L ⁻¹ dH ₂ O	5ml	2.05×10^{-9}
Cyanocobalamin (B ₁₂)	0.1g L ⁻¹ dH ₂ O	5ml	3.69×10^{-10}

(5) Sodium metasilicate	per 200 ml
Na ₂ SiO ₃ ·9H ₂ O	6 g
(Silicate can react with glass, ensure this stock solution is stored in plastic bottles).	

Medium	per litre
Stock solution 1	1.0 ml
Stock solution 2	1.0 ml
Stock solution 3 (Trace elements)	0.1 ml
Stock solution 4 (Vitamin mix)	1.0 ml
Stock solution 5	1.0 ml

Make up to 1 litre with filtered natural seawater. Adjust pH to **8.0** with 1M NaOH or 1M HCl prior to autoclaving. For agar, add 15 g per litre Bacteriological Agar. Autoclave at 15 psi for 15 minutes.

f/2 + Si Modified

Created on: 05 Nov 2019	CCAP (Culture Collection of Algae and Protozoa), SAMS Ltd, Scottish Marine Institute, Oban, Argyll, PA37 1QA, UK Tel: +44 (0)1631 559000 Fax: +44 (0)1631 559001 Email: ccap@sams.ac.uk Web: www.ccap.ac.uk	2 Pages
----------------------------	--	---------

This medium uses additional sodium metasilicate. The higher concentration of metasilicate has been observed to help some *Coscinodiscus* cultures to grow better and slow down size reduction.

To make this medium, omit stock 5 and instead add 0.1705 g of sodium metasilicate ($\text{Na}_2\text{SiO}_3 \cdot 9\text{H}_2\text{O}$) per litre of medium.

Reference

Guillard RRL & Ryther JH (1962) Studies of marine planktonic diatoms. I. *Cyclotella nana* Hustedt and *Detonula confervaceae* (Cleve) Gran. Can. J. Microbiol. 8: 229-239. – adapted for CCAP

Reviewed 26th January 2024