

## ASW 300 + Barley (Artificial seawater with salinity of 300 psu)

Stocks	per 1000 cm <sup>3</sup>
(1) NaCl	243.2 g
KCl	7.6 g
MgCl <sub>2</sub> .6H <sub>2</sub> O	54.4 g
MgSO <sub>4</sub> .7H <sub>2</sub> O	59.4 g
CaCl <sub>2</sub> .2H <sub>2</sub> O	1.3 g

### Medium

Add 20 ml of the autoclaved artificial seawater with 300 psu (1) and 4 autoclaved barley grains to a T-25 tissue culture flask (Corning 25100).

Used for CCAP 1959/1 *Pleurostomum flabellatum*