

CH (Chalkley's Medium)

Freshwater Protozoa

Stocks**per 100 ml**

(1) NaCl	2.0 g
(2) KCl	0.08 g
(3) CaCl ₂	0.12 g

Medium**per litre**

Stock solutions 1-3	5.0 ml each
---------------------	-------------

Make up to 1 litre with deionised water. Autoclave at 15 psi for 15 minutes.

Reference

Chalkley, H. W. (1930). STOCK CULTURES OF AMEBA. Science, 71(1843), 442-442. - adapted for CCAP

Reviewed: 17th January 2025