

EG (*Euglena gracilis* medium)

Freshwater algae and protozoa

Stock	per litre
-------	-----------

(1) CaCl ₂	1.0 g
-----------------------	-------

Medium	per litre
--------	-----------

Sodium acetate trihydrate	1.0 g
Lab-Lemco powder	1.0 g
Tryptone	2.0 g
Yeast extract	2.0 g
Stock solution (1)	10.0 ml

Add the above constituents and make up to 1 litre with deionised water. For agar, add 15 g per litre Bacteriological Agar. Autoclave at 15 psi for 15 minutes.

Reference

Schlosser UG (1982) Sammlung von algenkulturen. Berichte Der Deutschen Botanischen Gesellschaft 95: 181–276. – adapted for CCAP

Reviewed: 24th January 2024