

EG (*Euglena gracilis* medium)

Freshwater algae and protozoa

Stock per litre

(1) CaCl₂ 1.0 g

Medium per litre

Sodium acetate trihydrate 1.0 g
"Lab-Lemco" powder (Oxoid L29)* 1.0 g
Tryptone (Oxoid L42)* 2.0 g
Yeast extract (Oxoid L21)* 2.0 g
Stock solution (1) 10.0 ml

Add constituents above and make up to 1 litre with deionized water. For agar add 15 g per litre *Bacteriological Agar (Oxoid L11). Autoclave at 15 psi for 15 minutes.

Supply

* Unipath Ltd, Wade Road, Basingstoke, Hants RG24 0PW, UK

Reviewed: 6th August 2020