

 Culture Collection of algae & protozoa	QF 09 Media Preparation Recipe Template	Revision: 001 Approved: 30APR20
--	--	--

MW (Mineral Water)

Freshwater protozoa

Medium

There are many proprietary brands of bottled mineral water. At present, CCAP uses "**Volvic**". Dispense into suitable vessels and autoclave at 15 psi for 15 minutes.

Reviewed: 17th January 2025

Created on: 05 Nov 2019	CCAP (Culture Collection of Algae and Protozoa), SAMS Ltd, Scottish Marine Institute, Oban, Argyll, PA37 1QA, UK Tel: +44 (0)1631 559000 Fax: +44 (0)1631 559001 Email: ccap@sams.ac.uk Web: www.ccap.ac.uk	Page: 1 of 1
----------------------------	---	--------------