

 Culture Collection of algae & protozoa	QF 09 Media Preparation Recipe Template	Revision: 002 Approved: 29 MAY 26
--	--	--------------------------------------

NA (Nutrient Agar)

Medium

Nutrient agar

per litre

28.0 g

Thoroughly disperse the agar in deionised water in suitable vessels. Autoclave at 15 psi for 15 minutes.

Reference

Page, F.C. (1988) *A New Key to Freshwater and Soil Gymnamoebae*. Freshwater Biological Association.

Reviewed: 4th June 2026

Created on: 29 May 2026	CCAP (Culture Collection of Algae and Protozoa), SAMS Ltd, Scottish Association for Marine Science, Oban, Argyll, PA37 1QA Tel: +44 (0)1631 559000 Email: ccap@sams.ac.uk Web: www.ccap.ac.uk	Page 1 of 1
----------------------------	---	-------------