

PP (Proteose Peptone Medium)

Freshwater algae

Stocks (SES stocks)**per litre**

(1) K_2HPO_4	1.0 g
(2) $MgSO_4 \cdot 7H_2O$	1.0 g
(3) KNO_3	10.0 g

Medium**per litre**

Proteose peptone	1.0 g
Stock solutions 1 - 3	20.0 ml each

Make up to 1 litre with deionised water. For agar, add 15 g per litre Bacteriological Agar. Autoclave at 15 psi for 15 minutes.

Reviewed: 17th January 2025