

RPL:PJ

Freshwater protozoa

Medium

3:1 mixture

See separate recipes. Autoclave separately. Mix aseptically when cool.

RPL (Rye Grass-Prescott Liquid)

Medium

Dried milled rye grass *

Prescott's & James's Solution (PJ - see recipe below)

per litre

1.0 g

1.0 litre

Bring PJ to the boil then add rye grass. Continue to boil for 5 minutes. Allow to cool then restore volume to 1 litre with deionized water. Filter through GF/C paper and autoclave at 15 psi for 15 minutes.

Supply

* Please enquire

PJ (Prescott's & James's Solution)

Stocks

- (1) $\text{CaCl}_2 \cdot 2\text{H}_2\text{O}$
- KCl
- (2) K_2HPO_4
- (3) $\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$

per 100 ml

0.43 g

0.16 g

0.51 g

0.28 g

Medium

Stock solutions 1 - 3

per litre

1.0 ml each

Make up to 1 litre with deionized water. Autoclave at 15 psi for 15 minutes.