

QF 09 Media Preparation Recipe Template

Revision: 001 Approved: 30APR20

SES (Soil Extract with Added Salts)

Freshwater and terrestrial protozoa and marine algae

Stocks		per litre
	(1) K ₂ HPO ₄	1.0 g
	(2) MgSO ₄ .7H ₂ O	1.0 g
	(3) KNO₃	10.0 g
	(4) Soil Extract 2 (See SE2 recipe)	

Medium per litre

Stock solutions 1 - 3 20.0 ml each Stock solution 4 (SE2) 100.0 ml

Make up to 1 litre with deionised water. Autoclave at 15 psi for 15 minutes.

Reviewed: 17th January 2025